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# **⚠ Information for use - Please read carefully!**

## **HCCart**

### **Sodium bicarbonate for haemodialysis**

#### **Composition**

1 cartridge contains:  
Sodium bicarbonate 750 g

#### **Composition of the alkaline bicarbonate haemodialysis concentrate**

Unless otherwise prescribed, 750 g sodium bicarbonate for haemodialysis and water of a suitable grade (see notes under Composition of ready-to-use bicarbonate haemodialysis solution) yields approximately 7.5 l concentrate.

1000 ml of this concentrate contains:

Sodium bicarbonate	100 g
equivalent to: Na <sup>+</sup>	1190 mmol
HCO <sub>3</sub> <sup>-</sup>	1190 mmol

#### **Composition of ready-to-use bicarbonate haemodialysis solution**

Alkaline bicarbonate haemodialysis concentrate + water of a suitable grade + acidic bicarbonate haemodialysis concentrate yields after the ready-to-use bicarbonate haemodialysis solution.

For the electrolyte concentration and theoretical osmolarity of this ready-to-use solution, please refer to the package insert of the acidic bicarbonate haemodialysis concentrate used.

Freshly distilled water obtained under sterile conditions is the preferred medium for the dilution of bicarbonate haemodialysis concentrates.

However, purified water (Aqua purificata) may also be used if it meets the microbiological requirements of drinking water and complies with the notes below. If the water is demineralised, special attention must be paid to the possible presence of pyrogens.

**Notes:** Drinking water is not suitable for the preparation of bicarbonate haemodialysis solutions.

If purified water is used, particularly in the case of repeated haemodialysis, it is necessary to be aware of the possible presence of trace amounts of water treatment residues or chemical elements.

In particular, it is recommended that the aluminium, tin, mercury, zinc, fluoride, phosphate and sulphate levels in the water used are monitored; a maximum aluminium concentration of 10 mcg/l should not be exceeded. It is also desirable for the water used for haemodialysis not to contain any unbound chloride and ozone.

#### **Indications**

The indications below apply exclusively to the ready-to-use haemodialysis solution prepared from alkaline and acidic bicarbonate haemodialysis concentrates diluted as directed:

- Acute renal failure
- Chronic renal failure
- Hyperhydration
- Intoxications
- Compensation of the acid-base and electrolyte balances
- Adjustment of the blood/plasma/body temperature

#### **Contraindications**

- Hyperkalaemia with the use of a potassium-containing acidic bicarbonate haemodialysis concentrate;
- Hypokalaemia with the use of a potassium-free acidic bicarbonate haemodialysis concentrate;
- Intractable coagulopathy.

In patients with unstable circulation and/or unstable blood pressure, a different extracorporeal method of treatment may be indicated.

#### **Side effects**

During haemodialysis, side effects such as hypotension, nausea, vomiting and convulsions may occur.

#### **Dosage instruction**

Unless otherwise prescribed, mix alkaline and acidic bicarbonate haemodialysis concentrate with water of a suitable grade (see notes under Composition of ready-to-use bicarbonate haemodialysis solution) as directed (see Composition). Please refer to the package insert of the acidic bicarbonate haemodialysis concentrate used for details of the mixing ratio.

#### **Mode and duration of administration**

To produce an alkaline bicarbonate haemodialysis concentrate for extracorporeal bicarbonate haemodialysis or bicarbonate haemodiafiltration.

**Notes:** Use an alkaline bicarbonate haemodialysis concentrate only in combination with an acidic bicarbonate haemodialysis concentrate diluted as directed. Monitor the ready-to-use bicarbonate haemodialysis solution before use by a suitable gas analysis method. The duration of administration must be determined on an individual basis.

#### **Notes**

Do not use sodium bicarbonate after the expiry date printed on the container.

Open cartridge immediately before preparing the concentrate. Discard any surplus.

Dilute alkaline and acidic bicarbonate haemodialysis concentrates immediately before use.

Discard any surplus.

Do not use damaged cartridges.

Do not use if concentrate is not clear.

The concentration of the ready-to-use bicarbonate haemodialysis solution must be monitored carefully.

#### **Keep out of reach of children.**

When setting up or disconnecting, pay attention to the directions for the relevant haemodialysis monitor, as described in the enclosed instructions.

To be used only with devices that have been approved for the use of bicarbonate powder. In order not to endanger the life of the patient, dialysis treatment should (must) only be carried out by qualified personnel, i.e. individuals who can operate and are familiar with the instructions for the use of the relevant haemodialysis monitor.

Powder is non-pyrogenic.

Storage: Keep between +5°C and +40°C. Protect against frost.  
Cartridge: PP

Read instructions carefully before using HCCart.  
This 750g cartridge is sufficient for one dialysis treatment of 5-6 hours at a dialysate flow rate of 500 ml/min. and sodium bicarbonate concentration of 35 mEq/l.

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