


AsahiKASEI


# A.V. Fistula Needle Set



ASAHI KASEI KURARAY MEDICAL CO., LTD.



# A.V. Fistula Needle Set




## Representative product lineup

Code	Gauge (Color coding)	Needle length (mm - inch)	Back eye	Wing type	Tube length
AVF-BF2530B-15	15G	25 - 1"	With	Fixed	300mm
AVF-BR2530B-15	15G	25 - 1"	With	Rotatable	300mm
AVF-BF2530-15	15G	25 - 1"	Without	Fixed	300mm
AVF-BR2530-15	15G	25 - 1"	Without	Rotatable	300mm
AVF-BF2530B-16	16G	25 - 1"	With	Fixed	300mm
AVF-BR2530B-16	16G	25 - 1"	With	Rotatable	300mm
AVF-BF2530-16	16G	25 - 1"	Without	Fixed	300mm
AVF-BR2530-16	16G	25 - 1"	Without	Rotatable	300mm
AVF-BF3230B-16	16G	32 - 1 1/4"	With	Fixed	300mm
AVF-BR3230B-16	16G	32 - 1 1/4"	With	Rotatable	300mm
AVF-BF2530B-17	17G	25 - 1"	With	Fixed	300mm
AVF-BR2530B-17	17G	25 - 1"	With	Rotatable	300mm

\* For the availability of other lineup, please contact us.

## Coding system

**AVF-B F 25 30 BG -16**



<b>Gauge</b>	15G, 16G or 17G
<b>Sterilization</b>	"Blank": ETO, G: Gamma
<b>Back eye</b>	"Blank": without backeye, B: with backeye
<b>Tube length</b>	150mm or 300mm
<b>Needle length</b>	20mm, 25mm or 32mm
<b>Wing type</b>	F: Fixed, R: Rotatable

## AsahiKASEI

**Global Network:**  
**ASAHIKASEI KURARAY MEDICAL CO., LTD.**  
 1-105 Kanda Jinbocho, Chiyoda-ku, Tokyo 101-8101, Japan  
 TEL:+81-3-3296-3727 FAX:+81-3-3296-3722  
 URL:<http://www.asahikasei-kuraray-medical.com>

**ASAHIKASEI MEDICAL EUROPE GmbH**  
 Herriotstrasse 1, 60528 Frankfurt am Main, Germany  
 TEL:+49-69-66371 500 FAX:+49-69-666 5193

**ASAHIKASEI MEDICAL AMERICA INC.**  
 3570 Winchester Road, Suite 103, Memphis, Tennessee 38118, USA  
 TEL:+1-901-362-6105 FAX:+1-901-362-6180

**EU-Regulatory Representative:**  
**Tee & Howe IP Attorneys GmbH**  
 Ratiborer Str. 10, 31789 Hameln, Germany

**Manufacturer:**  
**Bain Medical Equipment(Guangzhou) Co., Ltd.**

CE 0197

Copyright ©2010 Asahi Kasei Kuraray Medical Co., Ltd. Printed in Japan No. 2010.3-6148