




PERITONEAL DIALYSIS CATHETERS AND KITS



Health Line International's PERITONEAL DIALYSIS CATHETERS are available in multiple designs, lengths, and configurations. Strategically placed side holes allow for increased flow and reduced omental occlusions, while large polyester cuffs secure tissue in-growth preventing infection. With the addition of complete insertion trays, Health Line International is providing patients and practitioners various options to meet their dialysis therapy needs.

PERITONEAL DIALYSIS CATHETERS AND KITS



Catalogue Number	Box Quantity	Carton Quantity	French/mm	Length (cm)	# of Cuffs	Description
Peritoneal Straight Catheters (Adult, Pediatric, Neonatal)						
PCS-42-1	5	20	14.7/4.9	42	1	Straight
PCS-42-2	5	20	14.7/4.9	42	2	Straight
PCS-47-2	5	20	14.7/4.9	47	2	Straight
PCS-31-1N	5	20	14.7/4.9	31	1	Straight Neonatal
PCS-32-2P	5	20	14.7/4.9	32	2	Straight Pediatric
PCS-37-2P	5	20	14.7/4.9	37	2	Straight Pediatric
Peritoneal Coiled Catheters (Adult and Pediatric)						
PCC-39.25-1P	5	20	14.7/4.9	39	1	Coiled Pediatric
PCC-56.5-2	5	20	14.7/4.9	57	2	Coiled
PCC-62-2	5	20	14.7/4.9	62	2	Coiled
Peritoneal Placement Kits						
PC-K	5	20				

Package Contains:

- 1 Peritoneal Silicone Catheter
- 1 Female luer
- 1 Male luer cap
- 1 Catheter clip

Peritoneal Catheter Placement Kit Contains:

- 1 Introducer needle, thin-wall, 18 gauge x 2 3/4" (7cm) length
- 1 "J" Double-ended guidewire, 0.038" O.D. x 27.5" (70cm) length
- 1 Dilator, 10 French (3.3mm) O.D. x 6" (15.2cm) length
- 1 Catheter clip
- 1 Spreadable sheath introducer, 16fr
- 1 Tunneling stylet
- 1 #11 Scalpel
- 6 Gauze sponges
- 1 Syringe, 10cc
- 1 Female luer
- 1 Male luer cap

For ordering information, please contact your country's authorized representative or distributor.

For a list of distributor's and country representatives please email customerservice@hlic.net.



Health Line International Corp.
803 N. 1250 W. #1
Centerville, Utah
USA | 84014
PH: 1- 877- 847- 4542
www.hlic.net

European Authorized Rep.
Obelis s.a.
Boulevard General Wahis 53
1030 Brussels, BELGIUM
Tel: +(32)2.732.59.54
Fax: +(32)2.732.60.03
CE: 0197