


UNI-PICC™

UNI-PICC™ The ONE Choice for Peripheral Access

The UNI-PICC™ Non-Ct PICC is an economical way of achieving the effectiveness of PICC access when power injection is not necessary. Made from polyurethane, the UNI-PICC™ is a great option for a variety of intravenous therapies.

Low profile hub allows for suture or sutureless securement



Simple grip luers are easy on gloves and maintain strength and crack resistance

Features and Benefits

Polyurethane Construction – reverse taper technology helps reduce post insertion bleeding while maintaining radiopacity and patient comfort.

Low Profile Hub – allows for better patient comfort and easy securement.

High Flow Rates – large gauge lumens allow for optimal flow for fluid delivery.

Standardized Trays and Components – Trays are organized with familiar industry components for clinician safety and user friendly procedural steps, including tearaway sheaths, echogenic safety needles, and premium guidewires.

Multiple Kit Configurations – UNI-PICC kits are available in Nursing, Basic, I.R. Long Wire, and Max Barrier tray configurations.

UNI-PICC Size Configurations

Single Lumen	Double Lumen
4fr	5fr
5fr	6fr

Standard Tray Components

	Quantity
UNI-PICC™ Non-CT PICC	1
Needle Free Valve	1 or 2
0.15 x 75cm Twisted Wire Stylet	1
.030 I.D. Adaptor W/Sideport	1
21GA X 7CM SAFETY SELD. ECHO.NEEDLE	1
PTFE Tearaway Introducer	1
Nitinol Straight Mandrel GW W/ Tungsten Tip	1
10cc Disposable Luer Lock Syr	1
Safety Scalpel #11	1
24" Paper Tape Measure	1
PICC Stylet TAG	1
PICC Patient Info Guide	1
PICC Chart Sticker	1
PICC IFU	1
Securement Device	1

For ordering information or for question about the Health Line products, please contact your Health Line representative or send us an email at info@hlic.net.



Health Line International Corp.
803 North 1250 West #1
Centerville, Utah 84014 USA
Phone: (801) 773-7798
Fax: (801) 820-8007
www.hlic.net

ISO13485
Certified Quality System